

Microlab 300

QUICK OPERATING GUIDE

Software Code	Designation	Working Reagent (WR) Preparation	Sample Blank Preparation	Working Reagent (WR) + Sample Preparation	Comment
ALBU	ALBUMIN	Ready to use	NA	WR (1000 µL) + Sample (10 µL) <i>Incubation time = 5 min</i>	
PASL	ALP DEA	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (20 µL) <i>Incubation time = 30 sec</i>	
ALSL	ALT / TGP	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (50 µL) <i>Incubation time = 30 sec</i>	
AMSL	AMYLASE	Ready to use	NA	WR (1000 µL) + Sample (10 µL) <i>Incubation time = no</i>	
ASSL	AST / TGO	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (50 µL) <i>Incubation time = 30 sec</i>	
BITD	BILIRUBIN DIRECT	Ready to use	R1 (1000 µL) + Sample (100 µL) <i>Incubation time = 5 min</i>	R1 (800 µL) + Sample (100 µL) <i>Incubation time = 5 min</i> Add R2 (200 µL) <i>Incubation time = 50 sec</i>	Cleaning after each sample: - 1x Cleaning cycle (Bleach solution 5%) - 1x rinse cycle (distilled water)
BITD	BILIRUBIN TOTAL	Ready to use	R1 (1000 µL) + Sample (50 µL) <i>Incubation time = 5 min</i>	R1 (800 µL) + Sample (50 µL) <i>Incubation time = 5 min</i> Add R2 (200 µL) <i>Incubation time = 5 min</i>	
CALA	CALCIUM	Ready to use	NA	WR (1000 µL) + Sample (20 µL) <i>Incubation time = no</i>	Be careful : This application requires a filter at 660 nm .
CHLO	CHLORIDE	Ready to use	NA	WR (1000 µL) + Sample (10 µL) <i>Incubation time = 5 min</i>	
CHSL	CHOLESTEROL	Ready to use	NA	WR (1000 µL) + Sample (10 µL) <i>Incubation time = 5 min</i>	
HDLL	HDL CHOLESTEROL	Ready to use	R1 (1200 µL) + Sample (10 µL) <i>Incubation time = 4 min</i>	R1 (900 µL) + Sample (10 µL) <i>Incubation time = 5 min</i> Add R2 (300 µL) <i>Incubation time = 4 min</i>	
LDLL	LDL CHOLESTEROL	Ready to use	R1 (1200 µL) + Sample (10 µL) <i>Incubation time = 4 min</i>	R1 (900 µL) + Sample (10 µL) <i>Incubation time = 5 min</i> Add R2 (300 µL) <i>Incubation time = 4 min</i>	
CMSL	CK-MB	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (40 µL) <i>Incubation time = no</i>	

Working Temperature 37°C. Before use, the working reagent should be brought to the working temperature.
For more details, consult the Microlab 300 User Manual.

Microlab 300

QUICK OPERATING GUIDE

Software Code	Designation	Working Reagent (WR) Preparation	Sample Blank Preparation	Working Reagent (WR) + Sample Preparation	Comment
CKSL	CK NAC	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (40 µL) Incubation time = no	
CRCO	CREATININE JAFFE	1 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (100 µL) Incubation time = no	
GISL	GAMMA-GT	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (75 µL) Incubation time = no	
GPSL	GLUCOSE	Ready to use	NA	R1 (1000 µL) + Sample (10 µL) Incubation time = 10 min	
LLSL	LDH-L	Ready to use	NA	R1 (1000 µL) + Sample (35 µL) Incubation time = 3 min Add R2 (250 µL) Incubation time = no	
MAGX	MAGNESIUM XYLIDYL	Ready to use	NA	WR (1000 µL) + Sample (10 µL) Incubation time = 5 min	
PRTU	MICROPROTEIN	Ready to use	NA	WR (1000 µL) + Sample (20 µL) Incubation time = 5 min	
PHOS	PHOSPHORUS	Ready to use	NA	WR (1000 µL) + Sample (10 µL) Incubation time = 5 min	
PROB	TOTAL PROTEIN	Ready to use	NA	WR (1000 µL) + Sample (10 µL) Incubation time = 10 min	
TGML	TRIGLYCERIDES	Ready to use	NA	WR (1000 µL) + Sample (10 µL) Incubation time = 10 min	
URSL	UREA	4 vol. R1 + 1 vol. R2	NA	WR (1000 µL) + Sample (10 µL) Incubation time = no	
AUML	URIC ACID	Ready to use	NA	WR (1000 µL) + Sample (25 µL) Incubation time = 5 min	

Working Temperature 37°C. Before use, the working reagent should be brought to the working temperature.
For more details, consult the Microlab 300 User Manual.