



Turbodyne™ SC

Turbidimetry- Nephelometry Analyzer

KEY FEATURES

- Dual technology Optics SECURES Turbidimetric/ Nephelometric measurements selectively.
- Lot specific assay Smart Card (SC) SECURES Master Calibration, Optimised Mixing and Accurate Reading of reactions.
- Touch screen Display SECURES user friendly operations with online reaction curve.
- QC software with Levey Jennings SECURES reliable reporting.
- Electronic Pipette SECURES Auto start of testing time.
- Optimix mixer SECURES automatic onboard mixing of reaction mixture.
- Curve Correction SECURES sample correlation.
- Unique cuvette design SECURES optimised onboard mixing.
- Disposable cuvettes SECURES no carryover effects, reduces labour and minimizes biohazard.

TECHNICAL SPECIFICATIONS

Principle	Turbodyne SC™ uses dual technology optics and is based on Turbidimetry- Nephelometry principle. A light beam with intensity (I_0) passes through the test cuvette with Ag-Ab reactions and falls on the photodetector placed at 180° and 135°. In turbidimetric assays reduction in the intensity of transmitted light (I_t) measured at 180° corresponds to the concentration of test analyte. In nephelometric assays, increase in the amount of scattered light (I_s) measured at 135° corresponds to the concentration of test analyte. The selection of test principle is through pre-calibrated SMART CARD.
Optical System	Laser LED (650nm) and Silicon detectors at 135° and 180°
Measuring Chamber	Temperature $37^\circ \pm 1^\circ$ Automixing
Electronic Pipette	For reagent auto start reading
Smart Card	Lot specific assay calibration with applications
QC	3 level with Levey Jennings plot

Display	4.3 inches Colour touch screen, online reaction curve
Printer	Inbuilt thermal (58nm)
Data Transfer	RS 232, USB
Data Storage	500 results
External Barcode (Optional)	RS 232, USB
Power Supply	100-240 V, AC 50/60Hz, 12 V Power adaptor
Battery	2200 mAH li-ion battery
Dimensions	290 mm(L) x 170mm (W) x 110mm (H)
Weight	Gross Wt -2 Kg, Net Wt. -1.5 kg
Cat. No.	826TD000000

LIST OF STANDARD COMPONENTS SUPPLIED ALONG WITH THE INSTRUMENT

Components	Qty	Components	Qty
1. Incubator	1 No.	6. Electronic pipette (20-200µl)	1 No.
2. Power Cable Adaptor	2 Nos.	7. Manual pipette (5-50µl)	1 No
3. RS232 Cable	1 No.	8. Manual pipette (100-1000µl)	1 No
4. User Manual	1 No.	9. Dust Cover	1 No.
5. Paper Roll	1 No.		