

TECHNICAL SPECIFICATIONS

Photometric Module

- 6 interference filters (EMF) 340, 405, 505, 546, 578 and 630 nm and 2 optional
- Band width (± 8 nm)
- Silicon Photo Diode detectors
- Photometric range from 0 to 2.5 O.D
- Halogen lamp of 12 V, 20 W

Measuring Module

- Dual mode (Flow cell & Cuvette)
- 30 μ L flow cell volume
- 10 mm square cuvette
- Variable temperature range through 25°C, 30°C, 37°C by means of Peltier controlled elements
- Minimum aspiration volume – 300 μ L

Analysis Mode

- 1-Point linear & non linear
- 2-Point linear & non linear
- 1-Point sample blank linear & non linear
- Rate linear & non linear
- Absorbance
- Immunoturbidimetry
- Bichromatic
- Kinetic with calibration
- Kinetic with factor

Calibration Curve

- Up to 9 calibration points
- Up to 3 replicates per point
- Axis: Linear and log

Quality Control

- Levey-Jennings control chart
- Westgard algorithm

Display

- Touch screen and popup key pad
- View Area 120 x 92 mm
- 25 test menu in each screen
- Touch pen provided for easy scratch proof operation of touch screen

Memory

- 200 defined programs
- 5000 Test Results
- Reagent Blank storage
- 18000 Control Results
- Built in Calibration History
- Help Function Key

Result Recall

- Collated Report by Date, ID
- Date and Test wise

Interface

- LIS & PC compatibility
- Port for external Key board, RS 232 Port

Power

- Power supply: 85 - 265 V AC, 50-60Hz, 90W

Dimension

- 362 x 364 x 195 (L x W x H)
- Weight: 6 kg

Distributed by:



ISO 9001:2008 & ISO 13485:2012

AGAPPE DIAGNOSTICS LTD.

Corporate off: / Factory: "Agappe Hills", Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.
TEL: + 91 484 2867000, FAX: +91 484 2867222, productcorp@agappe.in | www.agappe.com

YOUR BEST PARTNER IN DIAGNOSTICS



Company reserves the right to change any design and technical features of the product at any time. If needed.

ADU/BRQ/PMT/011/2015/R01



MISPAVIVA

The Synergy of Superior Technology & Quality

MISPAVIVA

Semi-automated Clinical Chemistry Analyser



One System

MISPA VIVA, true semi-automated clinical chemistry analyser is a symbol of superior technology and convenience, designed to deliver high accuracy, quality results with least turnaround time and low maintenance.

Three Superior Technology



Unique D^{MAX} Technology

In Kinetic Chemistry:

The unique feature that automatically recalculates abnormal sample results
Accurate results without rerun for high samples

Saves valuable time, reagent & money



PC compatibility

The convenient option to work through your computer

Adds to convenience in programming, sample testing, QC analysis, reporting etc

Conveniently sophisticated



Energy Matched Filter (EMF)

Ensures long filter life

Offers equal & stable intensity at all wavelengths for accurate estimation

Offers high degree of optimization

Perfection Predefined

Seven reasons for One choice

- 1 Embedded system with 5.7" TFT LCD color display with popup keypad**
Better menu clarity
Better visibility of reaction curve & data
Can switch off while at any menu or at any stage
User-friendly software to save time
Offers unparalleled convenience
- 2 Direct Access Key**
Direct access to test menu - ready for work at any time
Eliminates the need to memorize test menu number
Offers convenience and saves time
- 3 ERA Flow cell**
Ensures clear light passage - for accurate and reproducible results
Copper flow cell holder - quickly achieves stable temperature for optimal reaction
Specially designed 30 μ L flow cell - low sipping volume
Dual operating channel - Flow cell & cuvette
Offers unmatched accuracy economically
- 4 ALD Function**
ALD function to ensure extended life for the lamp
Significant power saving - low maintenance expense
- 5 Accuracy Enhancing Features**
High precision of 0.0001 resolution and excellent reproducibility
Aspiration air gap and Auto Empty Function after each sample
Offers enhanced Accuracy and rules out the carry over
- 6 5000 Test Results memory**
Effective retrieval system- by date, patient ID or test name
Enhanced storage and convenience
- 7 18000 QC result memory**
Designed with three levels of QC function
QC statistics displayed with Levy Jennings control chart with Westgard algorithm
High quality reporting

Ultimate Experience Guaranteed